

Circulating nucleosomes for Detection of Colorectal Cancer and High-risk Advanced Adenomas

Li-Chun Chang¹, Marielle Herzog², Nathalie Hardat², Dorian Pamart², Han-Mo Chiu¹

1.Department of Internal Medicine, National Taiwan University Hospital, Taipei 10002, Taiwan.

2.Belgian Volition SRL, 22 Rue Phocas Lejeune, Parc Scientifique Crealys, 5032 Isnes, Belgium


Background:

- Approximately **60%** of Fecal immunochemical tests (FITs) are **false-positive** and result in a lot of unnecessary colonoscopies.
- There is an unmet need for a simple test to supplement FIT to triage individuals at high risk who require colonoscopy referral.
- FIT is less effective in detecting **proximal** lesions even advanced neoplasms.
- We investigated the levels of circulating free nucleosomes containing different epigenetic modifications in patients referred for colonoscopy and the detection rate of proximal neoplasms.

Methods:

- 520 Patients enrolled including (i) CRC (n= 33), (ii) advanced adenoma (AA) (n=123, including 18 with AA>2cm); (iii) non-AA (n=168); (iv) non-neoplastic polyps (n=30); (v) colonoscopy negative controls (n=166).
- Plasma and stool samples were obtained prior to colonoscopy.
- Nucleosome levels were tested by using 7 different quantitative immunoassays (Nu.Q[®] assays; Belgian Volition SRL, Belgium) targeting H3.1-nucleosomes and different histone modifications.
- FIT was performed by using OC-SENSOR (Eiken Chemical Co., Ltd., Tokyo, Japan)

Circulating Nucleosome Levels in Combination with FIT Improves Detection of High-risk Advanced and Proximal Adenoma and Provides a Method to Reduce Unnecessary Colonoscopy




Author contact information:

1.National Taiwan University Hospital

Han-Mo Chiu: hanmochiu@ntu.edu.tw

2. Belgian Volition SRL

Marielle Herzog: M.Herzog@volition.com



Results:

- At a cut-off of 20µg/g, FIT had a sensitivity of 92.9% at 17.5% specificity for CRC. All the CRC patients and 83.3% of the high-risk AA (>2cm) patients were FIT positive.
- Combination of 2 Nu.Q[®] biomarkers: H3K27Me3 and H3R8Cit plus FIT in a logistic regression model showed improved sensitivity of 95% at 20.6% specificity, allowing the detection of all CRC patients and 94.3% of AA patients including all high-risk adenomas.

Diagnosis	Positive Test	
	FIT	Nu.Q [®] + FIT
CRC	33/33	33/33
Advanced Adenoma (AA)	112/123	116/123
<= 1cm	45/48	43/48
1-2 cm	52/57	55/57
> 2 cm (high risk AA)	15/18	18/18
Non-AA	141/168	149/167
Non-neoplastic polyps	28/30	25/30
Control	137/166	131/165

- FIT detected 46 out of 55 proximal adenomas whereas the combined model could detect 50 proximal adenomas including 3 proximal AA >2cm not being detected by FIT (p<0.05).
- Unnecessary colonoscopies are potentially reduced by 21.8% of the control and 23.3% of the non-neoplastic polyp subgroups compared to 17.5% and 6.7%, respectively with the FIT.