

# Advances in Liquid Biopsy for Glioblastoma Diagnosis and Monitoring through Nucleosome Epigenetic Modifications Tracking

Priscilla Van den Ackerveken<sup>1</sup>, Jonathan Decarpentrie<sup>2</sup>, Clotilde Hannart<sup>1</sup>, Nathalie Hardat<sup>1</sup>, Matteo Riva<sup>3,4</sup>, Nathalie Donis<sup>5</sup>, Fabienne George<sup>6</sup>, Lionel d'Hondt<sup>7</sup>, Jonathan Douxilfs<sup>2,5</sup>, Marielle Herzog<sup>1</sup>

1: Belgian Volition SRL, 22 Rue Phocas Lejeune, Parc Scientifique Crealys 5032 Isnes, Belgium, 2: Department of Pharmacy, Clinical Pharmacology and Toxicology Research Unit, Namur Research Institute for Life Sciences, Faculty of Medicine, University of Namur, Namur, Belgium. 3: Laboratory of Tumor Immunology and Immunotherapy, Department of Oncology, Catholic - University of Leuven, Leuven, Belgium. 4: Department of Neurosurgery, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium. 5: Qualiblood s.a., Research and Development Department, Namur, Belgium. 6: Biobank, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium. 7: Department of Oncology, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium.

## Glioblastoma: Challenges & Unmet Needs


- Glioblastoma (GBM) is the most aggressive primary brain tumor in adults.
- Standard management includes surgery, radiotherapy, and chemotherapy.
- High relapse rate and poor prognosis (median survival ~15 months).

➤ **Major gap: No circulating biomarkers validated for GBM.**

## Our Goal : Develop a liquid biopsy approach

- Develop a liquid biopsy approach for the diagnostic, the patient follow-up and early relapse detection of glioblastoma.
- **Investigating circulating nucleosomes and their epigenetic marks as new Biomarkers.**

## Methods :




Circulating H3.1 nucleosomes and their epigenetic marks (PTMs) levels were analyzed, using the Nu.Q® Immunoassays (Belgian Volition), in K2EDTA plasma from two independent cohorts including :

- **Cohort 1:** 67 plasma from GBM patients collected at diagnosis and 99 plasma from Healthy Donors
- **Cohort 2:** Multiple samples from 4 GBM patients collected longitudinally from diagnosis (pre-therapy; D-1) and throughout treatment, with a total of 21 samples assessed.

# Advances in Liquid Biopsy for Glioblastoma Diagnosis and Monitoring through Nucleosome Epigenetic Modifications Tracking


## Results at diagnosis : Epigenetic Profiling of Circulating Nucleosomes Enables GBM Detection



Box plot from Nu.Q® quantification (expressed as ng/mL) showing a significant increase of (A) the H3.1- (B) H3R8Cit-, (C) H3K4Me2-nucleosome levels, in GBM samples (n=66-67) compared to healthy samples (n=95 - 99). (D) Mutation in H3K27 is detectable in the plasma of GBM patient #3 (confirmed to be carrying this mutation by tissue biopsy) by the specific Nu.Q®H3K27M. p-values were determined by Mann-Whitney (\*p<0.05; \*\* p<0.01; \*\*\* p<0.001).

# Advances in Liquid Biopsy for Glioblastoma Diagnosis and Monitoring through Nucleosome Epigenetic Modifications Tracking

## Results « Monitoring » : Nu.Q®H3.1 Mirrors the Clinical Course of GBM Patients




Samples from four GBM patients were analyzed at the diagnosis (pre-therapy; D-1) and throughout treatment, with a total of 21 samples assessed. Nucleosome levels were subsequently analyzed using the Nu.Q® H3.1 immunoassay and H3.1-nucleosome concentration (ng/mL) have been reported in graphical form. Results were then compared with Magnetic Resonance Imaging (MRI) reports. RT (radiotherapy); CTX (Chemotherapy); FU (Follow-up); D (Days).



# Advances in Liquid Biopsy for Glioblastoma Diagnosis and Monitoring through Nucleosome Epigenetic Modifications Tracking

## Conclusions



- ✓ Histone mutation (H3K27M) in brain-GBM cells can be detected in blood samples using Nu.Q<sup>®</sup> Immunoassay (Nu.Q<sup>®</sup>H3K27M)
- ✓ High H3.1-nucleosome levels and the presence of H3K27M-nucleosomes in plasma from GBM patients suggest that nucleosomes cross the blood-brain barrier.
- ✓ The use of Nu.Q<sup>®</sup> PTMs tests on the GBM samples suggests that epigenetic markers may be useful for glioblastoma detection
- ✓ Circulating nucleosome levels mirrors the clinical course of GBM Patients

Altogether, our findings suggest that Nu.Q<sup>®</sup> immunoassays could serve as reliable, minimally invasive biomarkers for the detection and monitoring of disease progression in GBM patients

## Acknowledgements

This program benefited from a financial support of Wallonia in the frame of a BioWin's Health Cluster program. The illustrations were created using **BioRender** (<https://BioRender.com>).