

Priscilla Van den Ackerveken¹, Jonathan Decarpentrie², Clotilde Hannart¹, Nathalie Hardat¹, Matteo Riva^{3,4}, Nathalie Donis⁵, Fabienne George⁶, Lionel d'Hondt⁷, Jonathan Douxilfs^{2,5}, Marielle Herzog¹

1: Belgian Volition SRL, 22 Rue Phocas Lejeune, Parc Scientifique Crealys 5032 Isnes, Belgium, 2: Department of Pharmacy, Clinical Pharmacology and Toxicology Research Unit, Namur Research Institute for Life Sciences, Faculty of Medicine, University of Namur, Namur, Belgium. 3: Laboratory of Tumor Immunology and Immunotherapy, Department of Oncology, Catholic - University of Leuven, Leuven, Belgium. 4: Department of Neurosurgery, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium. 5: Qualiblood s.a., Research and Development Departement, Namur, Belgium. 6: Biobank, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium. 7: Department of Oncology, Université Catholique de Louvain, CHU UCL Namur, Yvoir, Belgium.

Glioblastoma: Challenges & Unmet Needs

- Glioblastoma (GBM) is the most aggressive primary brain tumor in adults.
- Standard management includes surgery, radiotherapy, and chemotherapy
- High relapse rate and poor prognosis (median survival ~15 months).
- > <u>Major gap</u>: No circulating biomarkers validated for GBM.

Our Goal : Develop a liquid biopsy approach

- Develop a liquid biopsy approach for the diagnostic, the patient follow-up and early relapse detection of glioblastoma
- Investigating circulating nucleosomes and their epigenetic marks as new Biomarkers.

Methods :



Circulating H3.1 nucleosomes and their epigenetic marks (PTMs) levels were analyzed, using the Nu.Q[®] Immunoassays (Belgian Volition), in K2EDTA plasma from two independent cohorts including :

- <u>Cohort 1:</u> 67 samples from GBM patients collected at diagnosis and 99 samples from Healthy Donors
- <u>Cohort 2</u>: Multiple samples from 4 GBM patients collected longitudinally from diagnosis (pre-therapy; D-1) and throughout treatment, with a total of 21 samples assessed.



Results at diagnosis : Epigenetic Profiling of Circulating Nucleosomes Enables GBM Detection

RespirERA

IHU Côte d'Azur

EUROPEAN LIQUID BIOPSY

SOCIETY



Box plot from Nu.Q[®] quantification (expressed as ng/mL) showing a significant increase of (**A**) the H3.1- (**B**) H3R8Cit-, (**C**) H3K4Me2-nucleosome levels, in GBM samples (n=66-67) compared to healthy samples (n=95 - 99). (**D**) Mutation in H3K27 is detectable in the plasma of GBM patient #3 (confirmed to be carrying this mutation by tissue biopsy) by the specific Nu.Q[®]H3K27M. p-values were determined by Mann–Whitney (*p<0.05; ** p<0.01; *** p<0.001).



Results Monitoring : Nu.Q®H3.1 Mirrors the Clinical Course of GBM Patients

RespirERA

IHU Côte d'Azur

EUROPEAN LIQUID BIOPSY

SOCIETY



Samples from four GBM patients were analyzed at the diagnosis (pre-therapy; D-1) and throughout treatment, with a total of 21 samples assessed. Nucleosome levels were subsequently analyzed using the Nu.Q[®] H3.1 immunoassay and H3.1-nucleosome concentration (ng/mL) have been reported in graphical form. Results were then compared with Magnetic Resonance Imaging (MRI) reports. *RT (radiotherapy); CTX (Chemotherapy); FU (Follow-up); D (Days).*





SMRE 2025

Conclusions



- Histone mutation (H3K27M) in brain-GBM cells can be detected in blood samples using Nu.Q[®] Immunoassay (Nu.Q[®]H3K27M)
- High H3.1-nucleosome levels and the presence of H3K27Mnucleosomes in plasma from GBM patients suggest that nucleosomes cross the blood-brain barrier.
- The use of Nu.Q[®] PTMs tests on the GBM samples suggests that epigenetic markers may be useful for glioblastoma detection
- Circulating nucleosome levels mirrors the clinical course of GBM Patients

Altogether, our findings suggest that Nu.Q[®] immunoassays could serve as reliable, minimally invasive biomarkers for the **detection** and **monitoring** of disease progression in GBM patients

Acknowledgements

This program benefited from a financial support of Wallonia in the frame of a BioWin's Health Cluster program. The illustrations were created using **BioRender** (https://BioRender.com).







